






**medicon**  
**Instrumente**




*Original*  
**GLASSMAN**



06.56.15 ●




06.56.24 ●  
06.56.20



06.58.15 ●

ORIGINAL GLASSMAN




06.58.24 ●  
06.58.20 20cm

GLASSMAN Tissue Forceps,  
jaws with non-crushing double-row  
serrations, tip 2 mm wide,


06.56.15 15 cm – 6"  
06.56.24 24 cm – 9 1/2"  
06.56.20 20cm

GLASSMAN Tissue Grasping Forceps,  
jaws with non-crushing double-row  
serrations,

06.58.15 15 cm – 6", jaws 4 mm wide  
06.58.24 24 cm – 6", jaws 8 mm wide  
06.58.20 20cm



50.03.15 ●  
50.03.20 ●



50.03.16 ●  
50.03.19 ●



**Medicon<sup>®</sup>**  
**Instruments**

## ORIGINAL GLASSMAN




50.11.16 ●  
50.11.20 ●  
50.11.23 ●


GLASSMAN-Allis Tissue and Intestinal  
Forceps, very delicate, jaws with  
non-crushing double-row serrations,  
50.03.15 15 cm – 6"  
50.03.20 20 cm – 8" (also useful as Common Duct  
Holding Forceps, see page 8)

GLASSMAN-ALLIS Tissue and  
Intestinal Forceps, modif. standard  
pattern with non-crushing double-row  
serrations,  
50.03.16 15 cm – 6"  
50.03.19 19 cm – 7 3/4"

GLASSMAN-BABCOCK Tissue  
Holding Forceps, narrow jaws with  
non-crushing double-row serrations,  
50.11.16 16 cm – 6 1/4"  
50.11.20 20 cm – 8"  
50.11.23 23 cm – 9"



50.30.24 ●




50.31.24 ●



**medicom<sup>®</sup>**  
**Instruments**

## ORIGINAL GLASSMAN

mit traumatischer GLASSMAN-  
Paralleleriefung.




**15.82.66**  
16,5 cm

**15.83.66**  
16,5 cm

**Bainbridge-Glassman**


GLASSMAN Gastro-Intestinal Clamps,  
jaws 65 mm long, with non-crushing  
double-row serrations,  
**50.30.24** 24 cm – 9 1/2", straight  
**50.31.24** 24 cm – 9 1/2", angular

GLASSMAN Gastro-Enterostomy  
Clamps, jaws 115 mm long, with  
non-crushing double-row serrations,  
**50.30.28** 28 cm – 11", straight  
**50.31.28** 28 cm – 11", angular



**50.30.28** ●

**50.31.28** ●





50.31.74

50.31.73


50.31.72

50.31.72 50 mm  
 50.31.73 60 mm  
 50.31.74 75 mm  
 24 cm  
**Glassman**



50.32.25 ●

## ORIGINAL GLASSMAN



50.33.25 ●

**50.31.74** GLASSMAN Anterior Resection Clamp,  
right-angled jaws 75 mm long, with  
non-crushing double-row serrations,  
24 cm – 9 1/2"

GLASSMAN Gastric Resection Clamp,  
jaws 105 mm long, with non-crushing  
double-row serrations

**50.32.25** 25 cm – 10", straight

**50.33.25** 25 cm – 10", angular

These clamps are also useful in wedge  
resection of the lung.